

HUGER 69 SERIES SMALL ANIMAL ENDOSCOPES

END-HGAGVE-SA-CONFIG*



Price: £8,400.00 - £9,261.00

Date Printed: 17/04/24

Category: Flexible Endoscopy > Diagnostic Imaging > Endoscopes > Video Endoscopes

PRODUCT DESCRIPTION

The 69 series of Huger endoscopes provide a full HD image when used alongside the HUGER 69 series Video processors.

Using this complete package, you will get an exceptional view of your procedure, with minimised noise and reduced halation. The full HD system gives you a viewing screen like no other. Capture and record images and videos at the click of a button.

Each endoscope is supplied in a hard shell transportation case and includes factory accessories.

Key Features

- Get a full-HD image resolution when used with the Huger 69 processors
- Tip Deflection: Up/Down:180 °/180 °, Left/Right:160 °/160 °
- Field of View: 145 °
- Keyboard included for easy patient data input
- Record still & video images to USB

5mmOD x 1050mmWL with 2.0mm channel Bronchoscope - Small animal HD video gastroscopy with 4 way angulation, commonly used for feline gastroscopy or small animal bronchoscopy

8mmOD x 1500mmWL x 2.2mm channel Videoscope - 4-way video gastroscope, commonly used for small animal upper and lower GI tract.

8mmOD x 1500mmWL x 2.8mm channel Videoscope - NEW 2.8mm channel, 4-way angulation. Commonly used for upper & lower GI tract and equine upper respiratory tract.

9.2mmOD x 1500mmWL x 2.8mm channel Videoscope - Commonly used for canine upper & lower GI tract and equine upper respiratory tract.

9.6mmOD x 1500mmWL x 2.8mm channel Videoscope - Commonly used for equine upper respiratory tract

12.9mmOD x 3300mmWL x 3.7mm channel Videoscope - Commonly used for equine gastroscopy

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	END-HGAGVE-SA-Config*
Procedure Type	Bronchoscopy, Cystoscopy, Respiratory Examination, Rhinoscopy
Working Length (mm)	1050mm
Outer Diameter (mm)	6.0mm

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights