

BURTONS K9 GASTROSCOPY TRAINING MODEL

END-TRM1



PRODUCT DESCRIPTION

The Burton's K9 Gastroscopy Training Model has been designed based on a 25-kilo dog, and contains a range of clinical disease states and pathology for veterinarians to investigate.

Made of synthetic material, the model is designed for flexible veterinary endoscopes up to 9.8mm in outer diameter.

Zippered entry enables the introduction/ removal of foreign bodies for retrieval training.

Ideal for Universities, Colleges and Veterinary Endoscopy Training centres, the model is supplied in a hard case with lubricant and can be used with/without the case for teaching purposes.

Anatomically the model is made up of:

- the oesophagus
- lower oesophageal sphincter
- fundus
- cardia
- body
- antrum
- pyloric sphincter
- duodenum

Disease states include:

- ulceration of the lesser curvature and thickened area
- petechial haemorrhage pyloric antrum
- peyers patches
- duodenal papilla
- bile stains
- nodular changes in duodenum

ADDITIONAL INFORMATION

SKU	END-TRM1
Brand	Burtons
Animal	N/A

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights