

PIEZON SCALER TIP TESTER CARD

118-205



PRODUCT DESCRIPTION

Check your scaler tips for wear with the Piezon Instrument Tester

With this scaler tip card, you will always be able to check any wear on your instruments, giving you peace of mind knowing your instruments are efficient and effective. This card will give you an indication of when to replace your instruments too. For example, you will have a 50% performance loss if your scaler tip is 2mm shorter than before.

To check your instruments, simply place the tip on the corresponding opening on the card.

ADDITIONAL INFORMATION

Animal	Mixed
SKU	118-205
Brand	Piezon
Autoclaveable	No

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights