

36INW X 24INH X 28.25IND DOUBLE DOOR BURTONS CAGE - CLEARANCE

170-028-2436-DD-CL



PRODUCT DESCRIPTION

Built from the highest quality materials, Burtons Lifetime Cages are guaranteed to last a lifetime!

36inW x 24inH x 28.25inD Burtons Double Door Cage

Hardwearing and easy to clean the Burtons 'Lifetime' range is the perfect option for everyday caging veterinary practices.

Burtons 'LIFETIME' stainless steel cages have been designed from the ground up to be strong and durable, easily standing up to the rigors of everyday practice life whilst still looking great. Supplied with sound deadening insulation as standard and constructed entirely in stainless steel for easy cleaning and increased hygiene (including the doors, hinges, latches and even screws!)

The cage body is made from the highest grade stainless steel and constructed for superior strength and elimination of dirt traps on seams. The thicker construction of the Lifetime cages reduces bends and un-shapping, that can cause the thunder noise when the patients move around.

All of our cages are available in standard sizes from stock with special made to measure sizes available on request. The completely modular, standardized design ensures not only an easy set up but allows your bank to grow with your practice. They are even compatible with any existing Burtons, Ken cage & Shorline cages that you may have.

Features

- Stainless steel, electro-polished removable door may be mounted either side
 - Covered rear corners on top and bottom eliminate accumulation of dirt and bacteria
 - Stainless steel sub frame makes the front of the cage extra strong
 - Radius bends on all sides provides a smooth, safe, easy to clean interior surface
 - Stainless steel self-latching door latch is designed for one handed opening and closing, with the ability to accommodate a padlock for greater security
 - Stainless steel doorframe and cross members are welded together for strength and rattle-less finish
- Distinctive floor design with fluid retaining lip as standard

- Stainless steel latch keeper latches door securely
- Heavy gauge stainless steel die formed barrel style hinges operate smoothly for the life of the cage
Latch connector bar is finger notched for smooth, easy opening
- Stainless steel hinge pins remove easily for cage maintenance
- Heavy 12 gauge stainless steel latch brackets are fully adjustable for smooth operation

S/N: 0000493074

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	170-028-2436-DD-cl
Brand	Burtons
Type	Cage
Complete Bank	No
Material	Stainless Steel
Mobile	No
Door Type	Reversible Mesh
Door Type	Double
Cage Width	36"
Cage Depth	28"
Cage Height	24"

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights