

**'C' MOUNT FOR 1/3IN CHIP SENSOR 0.2A 0.5X CAMERA ADAPTOR  
FOR MOTIC TRINOCULAR**

145-203

**PRODUCT DESCRIPTION****Complete your digital microscopy kit with the Motic C-mount adapter****Key Features**

C-Mount 1/3 Chip Sensor Adaptor - 0.2A 0.5x

Required to connect Motic Camera to Microscope

Compatible with Motic B1-223ASC Trinocular Microscope

**ADDITIONAL INFORMATION****SKU**

145-203

**Animal**

N/A

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights