

**EDAN VET SPO2 SENSOR CLIP - SMALL**

E-02.01.110507010

**PRODUCT DESCRIPTION****Edan Vet SpO2 Sensor Clip - Small**

SPO2 sensor clips for Edan Veterinary pulse oximeter sensors

**ADDITIONAL INFORMATION**

<b>Animal</b>	Mixed
<b>SKU</b>	E-02.01.110507010
<b>Type</b>	Sensor
<b>Monitor Fit</b>	Edan X12
<b>Accessory Parameter</b>	SpO <sup>2</sup>

**Good to know...**

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

