

## LABCOLD PHARMACY WALL MOUNTED/BENCHTOP GLASS DOOR REFRIGERATOR 66 LITRE

125-039-G-BUNDLE



Category: Spring Offers > Tour Products > Furniture & Lighting > Labcold Fridges > Fridges & Freezers > Medical > Pharmacy Refrigerators > Pharmacy Refrigerators > Refrigerators > Shop All

Price: £833.25 - £1,042.20  
Date Printed: 17/04/24

### PRODUCT DESCRIPTION

**Powered by the unique Intellicold controller, this professional 66 litre pharmacy refrigerator has been designed from the ground up for the safe storage of cold chain medicines**

Labcold IntelliCold Pharmacy Refrigerator now with NEW touchscreen controller. Designed from the ground up for the safe storage of cold chain medicines as advised by Public Health England.

With superior temperature uniformity, low energy consumption and extra functions, this refrigerator is ideal for storing your sensitive pharmaceuticals in perfect conditions, helping busy professionals monitor regulations and requirements.

This 66-litre wall mount or bench top fridge is compact and ideal for consult rooms with a minimum amount of space.

Featuring a 1 metre secondary temperature probe that can be placed anywhere inside the pharmacy fridge and a touchscreen controller which is completely flat so it can be thoroughly cleaned and disinfected. The controller displays alarms and records minimum and maximum temperatures for both temperature probes. It also records alarm states, fridge operation, door openings and temperature data for both probes at minute intervals.

You can only view and store this data on a computer by purchasing a specially encoded micro SD card. This data helps you with record keeping and regulatory inspections. Plus, if you have any issues, you can establish the exact duration and severity. You can find these in the related products below.

A **UKAS Temperature calibration** certificate gives you peace of mind that your temperatures are calibrated correctly. It can also be used if any health and safety issues arise around the correct storage of drugs and can be shown to an auditor during inspections.

Labcold are accredited to provide on-site UKAS calibrations for newly manufactured Labcold products, so you can order a UKAS temperature calibration as part of your purchase.

### Key Features:

- Fitted with the NEW IntelliCold® touch screen controller
- 24/7 automatic temperature data logging which is downloadable via optional micro SD card
- Controller and alarm battery back up
- Door open alarm
- Power failure alarm
- High temperature alarm
- Low temperature alarm
- Minimum/maximum temperature recording and display
- Off cycle auto defrost
- Fan assisted cooling
- Access port
- Environmentally friendly refrigerant
- Ammonia free
- Door lock with two keys
- Self-closing door
- Brackets for wall mounting
- Real time temperature display for two temperature probes
- Secondary (wandering) temperature probe
- Internal light
- Remote alarm contacts
- Optional micro SD card for transferring temperature data to your computer
- FREE five year parts and labour warranty

### Specifications:

- Configuration Wall mounted/Benchtop
- Capacity (Litres) 66 litres
- Height x Width x Depth 735mm x 450mm x 510mm
- Weight 35kg
- Shelves 3
- Temperature Range +2°C to +8°C
- Door Type Glass
- Lock Type Key lock
- Temperature Display
- High/Low Temperature Alarm
- Door Open Alarm
- Mains Failure Alarm
- Max/Min Temp. Recording
- Alarm battery backup
- Auto Defrost
- Access Port
- Remote Alarm Contacts
- Integral Temp. data logging

## PRODUCT OPTIONS & ADD ONS

### Fridge

1x Labcold Pharmacy Wall Mounted/Benchtop Glass Door Refrigerator 66 Litre

### Add-ons

1x Micro SD card for Labcold Refrigerators  
1x Labcold UKAS Temperature Calibration Service

## ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	125-039-G-BUNDLE
<b>Overall Product Dimensions (mm)</b>	735mm x 450mm x 510mm
<b>Colour/Style</b>	White
<b>Brand</b>	Labcold
<b>Alarm</b>	Yes
<b>Fridge Capacity</b>	66
<b>Door Type</b>	Glass

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights