

## CENTRIFUGE PRO-VET HE SUPPLIED WITH 5401 COMBINATION ROTOR

146-015



### PRODUCT DESCRIPTION

---

**Now the number 1 choice for most veterinary practices, the Burton's Pro - Vet HE is a compact, reliable and easy to use machine**

Supplied With 5401 Combination Rotor - Holds 24x Capillary Tube & 12x2ml Tubes

#### Key Features:

- Clarity & ease of use due to LED display of four memories
- Extra thick stainless steel bowl makes it easy to clean and rust free
- Strong and light thanks to its alloy and steel frame
- An exceptional quality brushless motor and industrial grade inverter ensure reliability.
- 4 programs of speed and time
- Quiet running, <60Db
- Tachometer speed certification
- Ambient mode allows for cooler running

#### Program / Speed / Time

1. Haematocrit (PCV) 13,000rpm 300Secs
2. Serum 1 6,000 rpm 600 Secs

3. Serum 2 5,000 rpm 300 Secs

4. Urine 2,500 rpm 300 Secs

**Safety Features:**

- A multi-point lock system, emergency release, spring strut and lock detection ensure complete lid safety.
- Eradicates user loading error using imbalance detection
- Over speed and overheat sensors

**Specifications**

- Dims 280x240x360mm
- Weight 12.8kg
- Power (watts) 140
- Rotor BRK5401 & BRK5424
- **2 Year Warranty Provided as standard**

|                          |                         |
|--------------------------|-------------------------|
| Rotor                    | BRK5401                 |
| Rotor Type               | 24 x capillary & 12x2ml |
| Tube size max            | 2 x 75mm & 11 x 40mm    |
| Minimum Rcf (G)          | 10                      |
| Maximum speed Rpm        | 13,000                  |
| Maximum Rcf (G)          | 15,800                  |
| Radius max cms sample    | 8.4                     |
| Tube angle               | 0 & 60                  |
| Acceleration time (secs) | 25                      |
| Deceleration time (secs) | 25                      |
| Autoclavable (frequency) | 121C (10)               |

**ADDITIONAL INFORMATION**

|  |                    |
|--|--------------------|
| <b>Animal</b>                          | N/A                |
| <b>SKU</b>                             | 146-015            |
| <b>Overall Product Dimensions (mm)</b> | 235H X 235W X 350D |
| <b>Colour/Style</b>                    | White              |
| <b>Brand</b>                           | Burtons            |

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights