

## CENTRIFUGE PRO-VET HE SUPPLIED WITH 5401 COMBINATION ROTOR

146-015



Category: Spring Sale > Tour Products > Diagnostic Imaging > Shop All Spring Offers > Laboratory > Centrifuge

Price: £1,105.10  
Date Printed: 07/07/22

### PRODUCT DESCRIPTION

**Now the number 1 choice for most veterinary practices, the Burton's Pro - Vet HE is a compact, reliable and easy to use machine**

Supplied With 5401 Combination Rotor - Holds 24x Capillary Tube & 12x2ml Tubes

#### Key Features:

- Clarity & ease of use due to LED display of four memories
- Extra thick stainless steel bowl makes it easy to clean and rust free
- Strong and light thanks to its alloy and steel frame
- An exceptional quality brushless motor and industrial grade inverter ensure reliability.
- 4 programs of speed and time
- Quiet running, <60Db
- Tachometer speed certification
- Ambient mode allows for cooler running

#### Program / Speed / Time

1. Haematocrit (PCV) 13,000rpm 300Secs

2. Serum 1 6,000 rpm 600 Secs
3. Serum 2 5,000 rpm 300 Secs
4. Urine 2,500 rpm 300 Secs

### Safety Features:

- A multi-point lock system, emergency release, spring strut and lock detection ensure complete lid safety.
- Eradicates user loading error using imbalance detection
- Over speed and overheat sensors

### Specifications

- Dims 280x240x360mm
- Weight 12.8kg
- Power (watts) 140
- Rotor BRK5401 & BRK5424
- **2 Year Warranty Provided as standard**

Rotor	BRK5401
Rotor Type	24 x capillary & 12x2ml
Tube size max	2 x 75mm & 11 x 40mm
Minimum Rcf (G)	10
Maximum speed Rpm	13,000
Maximum Rcf (G)	15,800
Radius max cms sample	8.4
Tube angle	0 & 60
Acceleration time (secs)	25
Deceleration time (secs)	25
Autoclavable (frequency)	121C (10)

### ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	146-015
<b>Overall Product Dimensions (mm)</b>	235H X 235W X 350D
<b>Colour/Style</b>	White
<b>Brand</b>	Burtons

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights