


LABCOLD PHARMACY ELECTRONIC GLASS DOOR REFRIGERATOR 340 LITRE

125-043-G



PRODUCT DESCRIPTION

Free standing IntelliCold® Pharmacy Refrigerator with NEW touch screen controller and key lock

A professional 340 litre pharmacy refrigerator powered by a unique IntelliCold® touch screen controller, designed from the ground up for the safe storage of cold chain medicines as advised by Public Health England in Protocol for Ordering, Storing and Handling Vaccines.

Key Features:

- Fitted with the NEW IntelliCold® touch screen controller
- 24/7 automatic temperature data logging which is downloadable via optional micro SD card
- Controller and alarm battery back up
- Door open alarm
- Power failure alarm
- High & Low temperature alarm
- Minimum/maximum temperature recording and display
- Off cycle auto defrost
- Fan assisted cooling
- Access port
- Environmentally friendly refrigerant
- Ammonia free
- Key lock with optional master key
- Self closing door
- Real time temperature display for two temperature probes
- Secondary (wandering) temperature probe
- Internal light
- Remote alarm contacts
- Optional micro SD card for transferring temperature data to your computer

- FREE five year parts and labour warranty
- Compliant with all leading professional guidelines

Specifications:

- Configuration Free Standing
- Capacity 340 Litres
- Height x Width x Depth 1510 x 600 x 700
- Weight 91kg
- Shelves 6
- Temperature Range +2°C to +8°C
- Glass Door Type
- Key lock
- Temperature Display
- High/Low Temperature Alarm
- Door Open Alarm
- Mains Failure Alarm
- Max/Min Temp. Recording
- Alarm battery backup
- Auto Defrost
- Access Port
- Remote Alarm Contacts
- Integral Temp. data logging

ADDITIONAL INFORMATION

Animal	N/A
SKU	125-043-G
Colour/Style	White
Brand	Labcold
Alarm	Yes
Fridge Capacity	340 litres
Door Type	Glass

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights