


## LABCOLD PHARMACY ELECTRONIC GLASS DOOR REFRIGERATOR 340 LITRE

125-043-G



### PRODUCT DESCRIPTION

---

## Free standing IntelliCold® Pharmacy Refrigerator with NEW touch screen controller and key lock

A professional 340 litre pharmacy refrigerator powered by a unique IntelliCold® touch screen controller, designed from the ground up for the safe storage of cold chain medicines as advised by Public Health England in Protocol for Ordering, Storing and Handling Vaccines.

#### Key Features:

- Fitted with the NEW IntelliCold® touch screen controller
- 24/7 automatic temperature data logging which is downloadable via optional micro SD card
- Controller and alarm battery back up
- Door open alarm
- Power failure alarm
- High & Low temperature alarm
- Minimum/maximum temperature recording and display
- Off cycle auto defrost
- Fan assisted cooling
- Access port
- Environmentally friendly refrigerant
- Ammonia free
- Key lock with optional master key
- Self closing door
- Real time temperature display for two temperature probes
- Secondary (wandering) temperature probe
- Internal light
- Remote alarm contacts
- Optional micro SD card for transferring temperature data to your computer

- FREE five year parts and labour warranty
- Compliant with all leading professional guidelines

### Specifications:

- Configuration Free Standing
- Capacity 340 Litres
- Height x Width x Depth 1510 x 600 x 700
- Weight 91kg
- Shelves 6
- Temperature Range +2°C to +8°C
- Glass Door Type
- Key lock
- Temperature Display
- High/Low Temperature Alarm
- Door Open Alarm
- Mains Failure Alarm
- Max/Min Temp. Recording
- Alarm battery backup
- Auto Defrost
- Access Port
- Remote Alarm Contacts
- Integral Temp. data logging

### ADDITIONAL INFORMATION

<b>SKU</b>	125-043-G
<b>Colour/Style</b>	White
<b>Brand</b>	Labcold
<b>Alarm</b>	Yes
<b>Fridge Capacity</b>	340 litres
<b>Door Type</b>	Glass

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights