

## U60 VET EDAN ULTRASOUND MACHINE

USEN-U60 VET-BUNDLE\*



Category: Winter Sale > Spring Offers > Imaging > Imaging > Ultrasound > Diagnostic Imaging > Ultrasound Systems > Shop All > Shop All

Price: £6,950.00 - £9,130.00  
Date Printed: 16/04/24

### PRODUCT DESCRIPTION

## Designed to deliver clear results, in a compact and comprehensive package, at a cost-effective price.

The U60 VET Ultrasound Machine provides clear imaging for confident diagnosis with a comprehensive set of applications for all patients, whilst keeping the unit compact and easy to use. The compact, lightweight design with integrated battery allows for easy manoeuvrability and use anywhere. Together, with a large 15" monitor, a backlit, easy-to-use control panel and dual transducer ports, the U60 Ultrasound Machine provides a highly efficient workflow. A high-quality colour image is the result of Tissue Harmonic imaging, spatial compound imaging, speckle reduction imaging and multi-beam technology to maintain high frame rates. The comprehensive software provides customisable user setting modification, quick image tools for optimisation, along with large-capacity data management.

This machine is supplied with 2 probes, Linear L12-5Q (5-12MHz) & Micro-convex (5-10MHz)

#### Key Features:

- Provides pan zoom, needle guide, dual live, triplex, PW, CW, HPRF and more
- Large 15" monitor for easy viewing
- Dual transducer ports for efficient workflow
- Multiple ports for data transfer
- Integrated battery and compact design for portability
- Tissue Harmonic Imaging for better detail display

- Spatial Compounding Imaging to improve contrast resolution
- Speckle Reduction Imaging to reduce image noise and clutter
- Multi-Beam Technology to maintain high frame rates.
- Backlit matrix keyboard
- Customisable user setting modifications.

## PRODUCT OPTIONS & ADD ONS

### Ultrasound Machine

1x U60 VET Edan Ultrasound Machine

### Additional Probes

1x Edan U60 Phased Array Probe 2.5MHz

### Additional Trolley

1x Edan U60 Luxury Mobile Trolley

## ADDITIONAL INFORMATION

Animal	N/A
SKU	USEN-U60 VET-BUNDLE*
Colour/Style	Light Grey
Ultrasound screen size	15 inch
Ultrasound colour	Colour

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights