

## EDAN AX3 VET ULTRASOUND MACHINE

USEN-AX3 VET-BUNDLE\*



### PRODUCT DESCRIPTION

**Designed to deliver uncompromised performance at a cost-effective price, The Acclarix AX3 VET Ultrasound Machine provides reliability alongside high specifications and additional features.**

The high spec resolution paired with TAI (Tissue Adaptive Imaging) results in superb, detailed imagery, particularly at depth, by continuously and automatically optimising the image. With a thoroughly thought-out design, the AX3 VET improves workflow efficiency whilst giving the clearest possible results. The Dual transducer ports will simplify workflow, whilst an extra screen above the keyboard provides a highly sensitive touch screen for easy setting modifications at your fingertips. All whilst still being ultra-light and compact for easy portability. The intelligent software provides dedicated pre-sets and measurement packages for veterinary applications, with a customisable user interface for personalised design, and multiple data management methods: DICOM, network, PDF, AVI, and BMP file export.

This machine is supplied with 2 probes, Linear L12-5Q (5-12MHz) & Micro-convex (5-10MHz)

#### Key Features:

- Superb, detailed resolution, particularly at depth
- Tissue Adaptive Imaging (TAI) continuously and automatically optimises imaging, allowing more focus on different kinds of animals.
- High-sensitive touch screen for daily operation
- Dedicated pre-sets and measurement packages for veterinary applications

- Customisable user interface for personalised design
- Multiple data management methods: DICOM, network, PDF, AVI, and BMP file export.
- Dual transducer design to improve diagnostic efficiency
- Unique dual-batteries design supports long continuous scanning, more than 2 hours
- Ergonomic control panel with sensitive touchscreen.
- Ultra-light and ultra-compact design for easy portability

## ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	USEN-AX3 VET-BUNDLE*
<b>Colour/Style</b>	Light Grey
<b>Ultrasound screen size</b>	15 inch
<b>Ultrasound colour</b>	Colour

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights