

## EDAN U60 PHASED ARRAY PROBE 2.5MHZ

USEN-E-02.01.214614



### PRODUCT DESCRIPTION

## Phased array probe for cardiac, compatible with the U60 Vet Machine

Edan U60 Phased Array Probe 2.5MHz

### ADDITIONAL INFORMATION

|                     |                     |
|---------------------|---------------------|
| <b>Animal</b>       | Mixed               |
| <b>SKU</b>          | USEN-E-02.01.214614 |
| <b>Colour/Style</b> | Light Grey          |

#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

