

## HUGER 69 FULL HD FELINE & CANINE FLEXIBLE ENDOSCOPY SOLUTION

END-HUGER69-SA-PROMO24



### PRODUCT DESCRIPTION

## The Huger 69 range of processors and scopes provide you with an exceptional view of your procedure

Minimised noise, and reduced halation, work alongside the full HD system to give you a viewing screen like no other. When used with the below 69 series range of scopes, the video centre provides a crisp image, with a dimmable light source and the ability to capture and record images at the click of a button.

#### Bundle Includes:

- Huger 69 video processor & light source
- Huger 69 Small Animal Videoscope 8mm OD x 1,500mm WL
- Huger 69 Ultra Slim Video Bronchoscope 3.2mm OD x 600mm WL
- Cytoscope 2.0mm, 187.5mm WL
- Askir Suction Unit

### ADDITIONAL INFORMATION

|                        |                        |
|------------------------|------------------------|
| <b>Animal</b>          | Small Animal           |
| <b>SKU</b>             | END-HUGER69-SA-PROMO24 |
| <b>High Definition</b> | Yes                    |

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights