

HUG-U-VAC POSITIONERS

HUG-BUNDLE*



PRODUCT DESCRIPTION

The HUG-U-VAC positioners conform gently to patient body contours, with minimal pressure, providing precise, comfortable, positioning

Designed for repeated use, the HUG-U-VAC conforms to your patient's body, with anatomically correct precision. Its material is engineered to reduce pressure sores, be MRI-safe, x-ray lucent, and radiolucent, whilst also conserving body heat.

To use, simply place your patient on a suitably sized positioner and use the hand vacuum pump to extract the air from the positioner, the material will mould to the patient's body. If you have a suction pump in your practice, this can be used to extract the air.

Positioners are sold separately unless purchased as part of the set. Vacuum pump sold separately / The complete set option includes 1 of each positioner, a vacuum hand pump and a duffel bag.

Key features

- Conforms to body contours.
- Pressure sore prevention
- Secure head and neck support
- Helps to conserve body heat.
- X-ray Lucent
- Radiolucent
- Easy to wipe down and disinfect.
- Can be used for pre-op preparation.
- Surgery-ventral, dorsal or lateral position.
- Post-op Recovery - prevent aspiration.
- Radiological Exams - precise positioning for accurate X-rays.
- Ideal for ultrasound and CT scans
- Aids emergency transportation

Positioner Options:

- HUG-H - Head Positioner
- HUG-S - Small Positioner - 16" x 24" - Puppies, Kittens
- HUG-M - Medium Positioner - 20" x 30" - Cats, small dogs
- HUG-L - Large Positioner - 24" x 36" - Medium & Large dogs, goats, sheep
- HUG-XL - Extra Large Positioner - 30" x 45" - Large dogs, goats, sheep
- HUG-J - Jumbo Positioner - 36" x 54" - large animals - alpacas, llamas, foals, mini horses

ADDITIONAL INFORMATION

SKU	HUG-BUNDLE*
Animal	Mixed
Overall Product Dimensions (mm)	16" x 24" - 20" x 30" - 24" x 36" - 30" x 45" - 36" x 54"
Colour/Style	Blue

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights