

LIVE CPD WEBINAR FOR YOUR WHOLE PRACTICE - CAPNOGRAPHY

CPD-WEB-CAPNO



PRODUCT DESCRIPTION

In veterinary education, high-level personalised CPD that is 100% relevant to your practice and the specific needs of your staff is hard to come by.

That is why, with Burton's Academy, we offer live practice webinar-styled education on the anaesthesia monitoring topics your team uses daily.

What will this webinar cover?

Capnography has become routine in many practices and has been a great step forward in patient monitoring, but there are some fundamental pitfalls.

In this live webinar, we discuss how to set up your capnograph, present solutions for hypo-and hypercapnia, and discuss how to troubleshoot certain waveforms. One common misconception is that carbon dioxide (CO₂) is the "bad guy", and we must get rid of it, but actually, too much or too little CO₂ in the blood creates pH imbalances which result in disturbances to normal homeostatic processes in the body.

So, when we monitor our patients for hypocapnia or hypercapnia, we're not only assessing the patient's ability to ventilate effectively, but we are using this information to ensure the levels of CO₂ are normal in the lungs, these levels are similar to the levels of CO₂ in the blood and maintain an optimal blood pH.

Capnography is so much more than waveforms...

What's included?

The session will run for approximately **1 hour**, with a Q&A session afterwards, and ongoing support offered via email and phone. If there are any particular topics you would like to run through, feel free to add them in the comments section of the booking link we will provide.

There is no limit on the number of people in your practice that can join the webinar, and a viewing link will be provided allowing those not on that shift to watch remotely. (It will not be available to download or watch on demand).

Written notes and CPD certificates are sent after the event.

Speaker:

Karen Heskin BVSc CertSAO MSc(VAA) MRCVS - Veterinary Clinical Consultant (Anaesthesia)

Dates & Times available: Once you have purchased your ticket online, we will send you a follow-up email with a link to choose a suitable date and to provide any remaining information that we require.

Sessions will be held on Tuesday afternoons which you can book up to 2 months in advance. A session lasts between 45-60 minutes, with 15 minutes of Q&A.

Please remember that attendance must be live, as the sessions will not be recorded. A Microsoft Teams link will be shared with those who wish to attend either in practice, or remotely.

ADDITIONAL INFORMATION

SKU	CPD-WEB-CAPNO
Animal	N/A

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights