

## SCRUB SINK 1 BAY WALL MOUNTED - CLEARANCE

115-245-CL



### PRODUCT DESCRIPTION

---

## Stainless Steel Scrub Sink 1 Bay Wall Mounted

Designed for ease of installation; back splash accommodates pipework so no wall chasing is required. Suitable for lever mixer taps (sold separately) Waste fitting supplied separately.

**FREE** Waste Fitting included.

Overall Dimensions:

**66 x 48 x 57 cm (W.D.H)**

### ADDITIONAL INFORMATION

**Animal**

Mixed

**SKU**

115-245-cl

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights