

## TONOPEN AVIA VET

117-655



### PRODUCT DESCRIPTION

**With over 40 breeds of dogs and cats possessing the predisposition to glaucoma, the Reichert Tono-Pen AVIA VET provides simplicity for monitoring this condition.**

Designed with the veterinary care professional in mind, the ergonomically designed Tono-Pen AVIA VET does not require daily calibration, has a long life, easy to replace the battery and can be used in any position.

#### Key Features

- **No Calibration** The Tono-Pen AVIA VET does not require daily calibration.
- **Accurate and Easy to Use** The Tono-Pen AVIA VET is designed to help eliminate operator error and functions in any position, making it ideal to use on all breeds of animals. Activated by the push of a button, the IOP measurement is displayed on two large, easy-to-read, LCD screens and can be used by right or left-handed users.
- **Intelligent** The Tono-Pen AVIA VET utilizes micro strain gauge technology and a 1.0mm transducer tip. The device displays individual readings and the average of 6 readings along with a statistical confidence indicator, ensuring accurate, repeatable, and reliable tonometry results.
- **Versatile** The Tono-Pen AVIA VET features the same proven accuracy as the clinically trusted Tono-Pen VET, but with a more ergonomic design, four times longer lasting battery, and larger LCD screens. The proven reliability of Reicherts technology allows operators to take IOP measurements confidently anytime, anywhere, and in any position.
- **Clinically Trusted** With the knowledge that IOP should be measured in all red eyes for which the cause is not immediately obvious, veterinarians worldwide have made the Tono-Pen brand a trusted part of their practice for over two decades. The proven reliability of the technology allows operators to take IOP measurements confidently anytime, anywhere, and in any position. The Tono-Pen AVIA VET will become the new industry standard.
- 20 years of clinical knowledge went into the development of the AVIA VET.
- Rapid, accurate and objective IOP measurement.
- Functions in any position.
- Portable and extremely lightweight.

- Large easy-to-read LCD display.
- Right or Left handed operation.
- Long life lithium-ion battery.
- Durable carrying case to protect your Tono-Pen during out of office visits

**Tono-Pen AVIA VET Starter Kit Includes:**

- Tono-Pen AVIA VET
- Carrying case
- 150 individually wrapped Ocu-Film + tip covers
- Tono-Pen AVIA VET Powercel Battery

**Specification:**

- Weight: 2.4oz. (71 g)
- Length: 6.25 in. (16 cm)
- Width: 0.75 in. (2 cm)
- Thickness: 1.75 in. (4.4 cm)
- Transducer: Micro Stain Gauge
- Measurement Range: 1 - 99 mmHg
- Power Source: Tono-Pen AVIA VET Powercel Battery Pack

**ADDITIONAL INFORMATION**

|               |              |
|---------------|--------------|
| <b>Animal</b> | Small Animal |
| <b>SKU</b>    | 117-655      |

**Good to know...**

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights