

BULB FOR WELCH ALLYN 2.5V OTOSCOPE HEAD

WA00200



Price: £15.00

Date Printed: 07/07/22

Category: Accessories > Diagnostic Bulbs

PRODUCT DESCRIPTION**Bulb For Welch Allyn 2.5V Otoscope Head**

fits models 20100, 21600, 38800, 26000, 41000, 235, 77900, 42700, 40510, 380, 385 & 395

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	WA00200
Brand	Welch Allyn
Diagnostic Instrument Type	2.8v Otoscope head

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

