

## BURTONS 3/8BSP O<sub>2</sub> PIPELINE LEAD 6 METERS - CLEARANCE

102-778-SP-CL



### PRODUCT DESCRIPTION

Burtons 3/8BSP O<sub>2</sub> Pipeline Lead 6 Meters

For Burtons Anaesthetic Machines  
6m Length

### ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	102-778-SP-cl
<b>Tubing Length</b>	6m
<b>MRI Compatible</b>	No

#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

