

**BULB FOR WELCH ALLYN 3.5V OTOSCOPE HEAD (WA03100)**

WA03100

**PRODUCT DESCRIPTION**

Bulb for Welch Allyn 3.5V otoscope head also fits models 20000, 25020, 21700, 73500, 20200, 20202, 28100, 27000, 27050, 26500, 41100, 43300, 25000, 25200, 20260

**ADDITIONAL INFORMATION**

<b>Animal</b>	Small Animal
<b>SKU</b>	WA03100
<b>Brand</b>	Welch Allyn
<b>Diagnostic Instrument Type</b>	3.5V otoscope

**Good to know...**

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights