

MICROBLUE BINOCULAR HEAD

147-200-2



PRODUCT DESCRIPTION

Interchangeable binocular head for MicroBlue Digital Microscope

- 45° inclined tubes
- Adjustable interpupillary distance between 48 and 75mm
- \pm 5 diopter adjustment on left tube
- Fully rotatable

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	147-200-2

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

