

**WHITE ECG ELECTRODE PADS - 3.6CM WITH POPPER X 50**

VSL-WS00S50

**PRODUCT DESCRIPTION**

**The Ambu WhiteSensor WS features a highly conductive solid gel with strong adhesion to ensure a good signal quality during short-term ECG monitoring applications.**

Thanks to the flexible foam backing material, the electrode is easy to apply.

**Key Features:**

- Highly conductive wet gel - immediate reliable recordings
- Comfortable foam backing - resistant to fluids
- 3 or 4 mm fitting connector - direct connection to cable, without adapters
- Radiolucent and MRI Conditional - can remain on the patient during X-ray and MRI procedures
- Convenient presentation of 5 pieces/strip - for fast application

**ADDITIONAL INFORMATION**

<b>Animal</b>	Small Animal
<b>SKU</b>	VSL-WS00S50
<b>Monitor Fit</b>	Vetronic Lightning, Vetronic Sapphire Plus, Vetronic Sentinel
<b>Accessory Parameter</b>	ECG

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights