

2.3MM OVAL SNARE MICRO - 20MM X 2350MM L (END-1102) - CLEARANCE

END-1102-CL



PRODUCT DESCRIPTION

Flexible Oval Snare 20mm - 2.3mm D x 2350mm L

Highly reliable and easy to maintain, our flexible endoscopic instruments offer consistent performance and are CE and FDA approved.

- Constructed of special corrosion resistant steel alloys
- Autoclavable handles
- Colour coded and marked for identification

ADDITIONAL INFORMATION

| | |
|-----------------------|---|
| Animal | Mixed |
| SKU | END-1102-CL |
| Procedure Type | Bronchoscopy, Cystoscopy, Respiratory Examination, Rhinoscopy |
| Channel Size | 2.3mm |
| Working Length | 2350mm |

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights