

2.3MM OVAL HEXAGONAL SNARE - 2350MM L - CLEARANCE

END-1106-CL



PRODUCT DESCRIPTION

Highly reliable and easy to maintain, our flexible endoscopic instruments offer consistent performance and are CE and FDA approved

- Constructed of special corrosion resistant steel alloys
- Autoclavable handles
- Colour coded and marked for identification

ADDITIONAL INFORMATION

Animal	Mixed
SKU	END-1106-cl
Procedure Type	Bronchoscopy, Cystoscopy, Respiratory Examination, Rhinoscopy
Channel Size	2.3mm
Working Length	2350mm

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights