

3.5V OPHTHALMOSCOPE HEAD

WA11710



Category: Diagnostic Essentials

Price: £176.50
Date Printed: 29/03/24

PRODUCT DESCRIPTION

Standard Welch Allyn 3.5V Halogen Ophthalmoscope Head

Key Features

- 6 apertures for all general and specialist use
- 28 lenses -25 to +40 diopters for better resolution
- Rubber brow rest prevents scratching of eyeglasses
- Illuminated lens dial clearly identifies the setting

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	WA11710
Size	X-Small
Colour/Style	Silver
Brand	Welch Allyn
Digital	N/A
Diagnostic Head Type	Ophthalmoscope

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights