

## COOLVIEW CLED23 T1 CEILING MOUNTED EXAMINATION LIGHT - CLEARANCE

CLED23T1C-CL-2



### PRODUCT DESCRIPTION

---

**Coolview CLED23 delivers the highest lighting performance and colour corrected light, reducing eye strain making the focus of the light vivid and tissues easily identified.**

#### Highest Light Quality

CLED23TM emits 60,000 lux of intense bright light at 0.5m, providing a pool of clear light during procedures. The light intensity can be altered through 4 easily selectable settings providing a suitable level of intensity for all procedures.

#### Large Illuminated Area

A large illuminated area of 210mm diameter makes it easy to examine large areas of tissue and makes it easier to target the light on the patient.

#### High Colour Rendition

CLED23 reproduces all of the visible colours naturally and accurately. The high Ra of > 93 describes the exceptionally high quality of white light whilst the high R9 of > 95 measurement describes how it accurately reproduces visible red.

#### Colour Control Settings

CLED23TM allows the user to adjust the colour temperature of the light source. This enhances the clinician's ability to recognise and distinguish tissue conditions. The warm white (3,500°K) setting is ideal for dermatological examinations and the daylight white setting (4,500°K) is ideal for the recognition of blood vessels.

## Environmentally Friendly

16 powerful and long-lasting LEDs ( $\leq 30,000$  hours) provide a maintenance-free light source whilst using minimal power consumption. The LEDs do not contain Mercury, Lead, Cadmium, CFCs, PoPs, VOCs or Halogens.

## Full Rotation of Light Head

The Coolview CLED23™ light head boasts 360° rotation. This allows the user to position the light beam to suit any procedure, without any restrictive stops or twisted cables. This also means the examination light can be used to provide ambient room lighting by reflecting light from the ceiling.

## Full Rotation at Ceiling Pole

Ceiling mounted models (CLED23T1C) also have unrestricted 360° rotation at the ceiling pole providing even more positioning flexibility.

## Long Arm Reach & Large Angular Movements

The exceptionally long, 1200mm arm reach of the CLED23 ensures the powerful beam can be positioned all around the patient. This allows for full examination of the patient without the need for patient repositioning during an examination.

## Spring-Balanced Positioning

Innovative spring-balanced arm system ensures the light is positioned effortlessly, remains fixed in place and is maintenance-free.

## Fully Enclosed Light Head

The light head is sealed to prevent the ingress of dust and moisture to IP43

## Intuitive Control Panel

The easy to use and ergonomic control pod allows all lighting features to be adjusted with a simple press of the button.

## Specification

- Light Intensity@ 0.5m : 60,000 Lux
- Illuminance Field : 190mm
- LED Life: 30,000 hours L70 B50
- Colour Temperature: 3,500°K & 4,500°K
- Base of arm rotation: 240°
- Head rotation: 360°
- Type of light source: LED
- Average bulb life (Hours): 30,000 L70 B50
- Light intensity (Lux): 60,000 @ 0.5m
- Field size: 190mm
- Colour Rendition Index (Ra): >93
- Colour Rendition Index (R9): >95
- Dimming Controls: 4 settings 10%-100%
- Number of LEDs: 16
- Power Consumption: True Power 21W
- Apparent Power: 20 VA
- Temperature rise : < 1°C
- IP Rating: 43

## Mechanical Detail

- Arm length: 1200mm
- Balance of articulated joints: Directional Load Compensation (DLC)
- Controls: Touch panel on rear of lamp head

S/N: 32577/19

## ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	CLED23T1C-cl-2
<b>Colour/Style</b>	White
<b>Brand</b>	Brandon
<b>Autoclavable Handle</b>	No
<b>Battery Backup</b>	Optional
<b>Colour Control</b>	Yes
<b>Colour Temperature (k)</b>	3,500ÅK - 4,500ÅK (2 Steps)
<b>Focusable</b>	No
<b>Head Type</b>	Single
<b>LED Life (hours)</b>	30,000 hours
<b>Light Source</b>	LED
<b>Lux @1 meter</b>	30,000
<b>Proceduredure</b>	Dental, Exam, Minor-Operation

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights