

## ELITE RETINOSCOPE STREAK

WA18245



### PRODUCT DESCRIPTION

---

**The Welch Allyn Elite Retinoscope makes refractions easier than ever by producing one of the brightest, sharpest, and cleanest streaks available.**

Enhanced light output makes refractions easier than ever, by producing one of the sharpest, brightest, and cleanest streaks available on the market. Discover unparalleled performance in a compact design with convenient, user friendly features.

#### Key Features

- Provides 2.8X the brightness of our standard retinoscopes
- Superior streak quality enables easier and quicker refractions even in patients with small pupils
- Crossed linear polarizing filter eliminates bothersome glare from trial lenses, helping you perform a more precise exam
- External focusing sleeve improves comfort and maneuverability
- Sealed optics keep the lenses cleaner for more effective exams
- Fits all Welch Allyn 3.5 V power sources (Handle Not Included)

#### ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	WA18245
<b>Overall Product Dimensions (mm)</b>	17.2cmW x 4.7cmH x 4.7cmD
<b>Generic Size</b>	X-Small

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights