

GALAXY VISITOR CHAIR

225-706-CONFIG



PRODUCT DESCRIPTION

Stylish yet practical waiting room seating solutions with the ergonomically designed Galaxy Visitor Chairs

Much consideration has been given to the design of the Galaxy Visitor Chairs, providing a functional, economical seating solution without compromising on quality, style or practicality.

Visitors and staff alike will appreciate the comfortable deep foam seats with rounded corners. The plastic feet prevent the chair from involuntarily moving across smooth practice floors, keeping them well-protected from scratches and scrapes.

The Galaxy Visitor Chairs are fully upholstered in anti-bacterial vinyl, providing durable, comfortable and practical seating that can be easily yet highly effectively cleaned, reducing cross-contamination risks within your practice.

Available with arms or without, the Galaxy Visitor Chairs are supplied fully assembled ready for immediate use in your practice. With 13 colour options to choose from, you can be certain to find the ideal seat to match your practice's design. Black and navy blue colour versions of the seats are available for direct order through our website, but all colours are available for order – simply contact Burtons and a member of our team will be happy to help you!

Chair linking kits are available for chairs with both arms and without, and individual chairs can also be safely stacked to a maximum of five high, offering a valuable space-saving options where needed.

A **5-year warranty** covers failure of all mechanical parts during normal use.

Key Features:

- Deep foam upholstery
- Rounded seat corners
- Safely stackable (maximum 5 high)
- Plastic feet to ensure stability and floor protection
- Anti-bacterial material
- Supplied fully assembled
- Chair linking kit available for both arm and no-arm seat versions

- 5-year warranty

Specifications:

- Dimensions: 82cm (height) x 54cm (width) x 52cm (depth)
- Seat Height: 43.5cm
- Maximum user weight: 110kg (17 stone)

ADDITIONAL INFORMATION

Animal	N/A
SKU	225-706-config
Brand	Sunflower Medical
Number Of Seats	1
Armrest	No

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights