

## SET OF 5 REUSABLE EAR SPECULA: 2MM 3 MM 4MM 5MM & 9MM (SHORT)

WA22100



### PRODUCT DESCRIPTION

## Reusable Ear Specula Set 2mm, 3mm, 4mm, 5mm & 9mm (Short)

- Material is lightweight, yet very durable
- May be autoclaved, boiled or cleaned with standard germicides
- Colour: Green
- Material: Polypropylene

### ADDITIONAL INFORMATION

|               |              |
|---------------|--------------|
| <b>Animal</b> | Small Animal |
| <b>SKU</b>    | WA22100      |

#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

