

VETARIO T30M INTENSIVE CARE UNIT

114-145



PRODUCT DESCRIPTION

Multi-purpose small animal Intensive Care Unit, providing the warmth and isolation essential for critical care needs for a variety of young, sick, injured or weak animals and birds.

With an internal floor area of around 24cm square, the unit is ideally suited to veterinary practices where portability and ease of cleaning are key benefits.

The T30 Advance is portable and stackable. It operates at 12v for safety and can be powered from the mains with the adaptor (supplied) or from a car cigarette lighter (adaptor available separately, part number: 15.45). The clear door is made of tough, UV resistant polycarbonate that can be closed and latched with one hand.

The digital control system shows temperature in either °C or °F and includes dual adjustable alarms to warn of variations of temperature inside and outside the cabinet beyond user-set thresholds. Simply alter the temperature in 0.1°C increments by using the digital controls to adjust between room temperature and 38°C. The digital controller guards against accidental temperature adjustment. Temperature stability is maintained by a small fan which runs at low speed to minimise noise and disturbance. The adjustable vent allows plentiful ventilation in warm rooms or at low set temperatures which is important for larger animals.

Cleaning the T30 is easy: the hinged door lifts off and the working parts and fascia simply unclip and slide out – no tools needed – leaving the cabinet empty for deep cleaning. To reduce the risk of infection Polygiene Biomaster™ Antimicrobial Protection is embedded within the incubator cabinet plastic during manufacture which prevents the growth of harmful bacteria.

Key features:

- Portable – integral carry handle
- Stackable
- Clear door for easy observation
- Easy to clean
- Quick warm-up

- Quiet, low speed fan
- Digital control system reading °C or °F
- Double-walled cabinet for heat and sound insulation
- Temperature range adjustable from room temperature to 38°C in 0.1°C increments.
- Dual alarms to warn of internal and external temperature fluctuations.
- Protected by Polygiene Biomaster antimicrobial additives
- Safe (12v operation)
- Simple to disassemble for deep cleaning
- Servicing options outside of guarantee period
- Made in the UK

Protected by Polygiene Biomaster antimicrobial additives - to reduce the risk of infection Brinsea incorporate Polygiene Biomaster Antimicrobial Protection into our Vetario T30 cabinets during manufacture. Polygiene Biomaster Protection is embedded within the plastic itself which prevents the growth of harmful bacteria helping to provide the optimum environment for young animals. Proven to quickly reduce bacterial growth by up to 99.99%

Specifications:

- Overall height - 297mm
- Overall width - 270mm
- Overall depth - 275mm
- Floor area - 240mm x 240mm
- Effective volume - 9.5ltr
- Power consumption - Max 36, Typical 18 watts
- Power supply 115v, 60Hz or 230v 50Hz as ordered
- Weight - 2.1kg

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	114-145

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights