

## ECG ELECTRODE PADS - 3CM WITH BANANA PLUG SOCKET -50 PACK - CLEARANCE

VSL-M00A50-CL



### PRODUCT DESCRIPTION

---

**The small Ambu BlueSensor M features wet gel, offset connector and suitable adhesion to ensure stable signal quality during short-term applications.**

Thanks to the occlusive backing material, the electrode is resistant to liquids.

Connection: 4mm fitting

**Short expiry date**

#### Key Features

- Offset connector mitigates the impact on signal quality from external forces like pulling of lead wires and pressure on electrodes
- Offset connector allows the user to connect and disconnect the lead wire without applying pressure to the patient or the sensor area, which reduces the risk of gel spreading and adhesive detachment from skin.
- Wet gel optimises contact between patient and electrode
- Ag/AgCl sensor ensures stable signal quality

### ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	VSL-M00A50-cl
<b>Type</b>	Other
<b>Monitor Fit</b>	Vetronic Lightning, Vetronic Sapphire Plus

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights