

EDAN HANDHELD PULSE OXIMETER CHARGING BASE

E-VEH100B-POD



PRODUCT DESCRIPTION

Edan Pulse Oximeter Battery Charger Stand Kit (UK)

Kit Includes:

- 1 x Ni-MH Battery
- 1 x Charger Stand with Power Cord

1 year warranty

ADDITIONAL INFORMATION

| | |
|---------------------------|-----------------------|
| Animal | Mixed |
| SKU | E-VEH100B-POD |
| Monitor Parameters | SpO ₂ |
| Screen Type | Digital Display |
| Power Source | Battery, Rechargeable |

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights