

**HEINIGER BLADE, SIZE 50**

HNG-707-905-A

**PRODUCT DESCRIPTION****Heiniger Snap-On #50 / 0,2 mm Blade**

For veterinary use only. Interchangeable blade set for Heiniger Opal.

Heiniger clipper heads are compatible with clippers from Andis, Moser, Oster and Wahl.

0.2mm

**ADDITIONAL INFORMATION**

<b>Animal</b>	Small Animal
<b>SKU</b>	HNG-707-905-A
<b>Clipper Speed(s)</b>	N/A
<b>Number of Clipper Speeds</b>	1
<b>Compatible with other brand blades ?</b>	Yes
<b>Replaceable Battery</b>	N/A

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights