

ADAPTABLE ANAESTHETIC MACHINE O2 ONLY WITH ALARM - CLEARANCE

102-062-A-CL



PRODUCT DESCRIPTION

AAA (Adaptable Anaesthetic Apparatus) a versatile and adaptable frame design which can be used to provide the basis of a variety of anaesthetic machines. The all stainless steel frame is fitted with two shelves approximately: 290mm x 450mm and is mounted on 75mm anti-static castors. The frame accommodates two F size cylinders which are secured within hinged cylinder clamps on the rear of the frame.

Oxygen only version is supplied with two bullnose single stage regulators fitted with Schrader quick release valves to enable rapid cylinder change. The unit can be used in conjunction with cylinder gas or connected directly to a piped gas supply without the need of any further modification.

The back bar is wide enough to accommodate two vaporisers (vaporisers sold separately). The O2 flush features a combined corner unit with emergency oxygen flush, non-return valve and safety pressure relief valve. Gas flow is controlled by a high accuracy rotameter, ideally suited for low flow anaesthesia.

Key Features:

- Back bar accommodates two vaporisers
- Holds two 'F' Size cylinders
- Rapid Cylinder change
- High accuracy rotameter 10lpm O2/121pm N2O
- Emergency O2 Flush
- Safety pressure relief valve & gas non return valve
- Piped gas ready
- Vaporiser sold separately
- Stainless Steel Frame
- Overall Size 46 x 57 x 138cm WxDxH

Technical Specifications	
Overall Size:	57cmL x 46cmW x 138cmH
Patient Circuit Connector Outlet	BS22mm Male
Cylinder Capacity	2 x 'F' Cylinders
Flow Control	4O2 0.1 10lpm and N2O 0.1 to 121pm

ADDITIONAL INFORMATION

Animal	Mixed
SKU	102-062-A-CL
Gas Type	O ²
Gas Supply	'F' Size Cylinder(s), Piped Gas
Flowmeter Type	Graduated Rotameter
Concentrator	Optional
Optional O2 Alarm	Yes
Optional N2O	Yes
Emergency O2 Flush	Yes

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights