

**LIGHT GUIDE CABLE 1.8M L X 3.5MM OD WITH ADAPTOR AE6 & AL7**

END-402-11010-CL

**PRODUCT DESCRIPTION**

**The light guide is 1.8m Long x 3.5mm OD and is suitable for Halogen light sources and 180W Xenon light source.**

AE (adaptor endoscope side)  
AL (adaptor light-source side)

AE 6 (MGB / Wolf)  
AL 7 (MGB / Wolf)

**S/N: 447682****ADDITIONAL INFORMATION**

<b>SKU</b>	END-402-11010-cl
<b>Brand</b>	MGB
<b>Animal</b>	Mixed
<b>Colour/Style</b>	Grey

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights