

EDAN HAND HELD PULSE OXIMETER PROTECTIVE COVER

E-VEH100B-COVER



PRODUCT DESCRIPTION

Edan Hand Held Pulse Oximeter Protective Cover. Change when broken. Compatible with Edan VE-H100B.

Blue Edan VE-H100B Hand Held Pulse Oximeter Protective Cover

The Edan VE-H100B protective cover is a soft, durable rubber case designed to shield your Edan H100B Handheld Pulse Oximeter from scratches and bumps. Its anti-slip design allows a secure fit and the soft silicone padding ensures great protection. It is easy to attach and replace, and offers reliable protection for everyday use. The unit can still be charged on the charging station with the cover on.

ADDITIONAL INFORMATION

Animal	N/A
SKU	E-VEH100B-COVER
Accessory Parameter	Other
Monitor Compatibility	Edan VE-H100B
Change Frequency	When Broken

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights