

BURTONS 5% TEC 3 ISOFLURANE VAPORISER - BODY ONLY - CLEARANCE

600-001-CL



PRODUCT DESCRIPTION

Burtons manufactured Vaportec Tec3 Isoflurane Vaporiser

A temperature-compensating, out-of-circuit key filled vaporiser in a much lighter design, familiar to most veterinary practitioners. Calibrated specifically for the specified anaesthetic drug using a laser refractometer.

Designed to be used in continuous flow "out of circuit" systems and give a per-set volume to volume concentration of vapour; to be self-compensating over the average working temperature range, and non-critical of input gas flow rates.

Fitting Type

Choose from the standard taper mounted fitting (cage mount) or the quick release option; allowing you to easily swap the vaporiser to different anaesthetic machines throughout the practice. Additional quick-release compatibility blocks will be required for each anaesthesia machine. (Available separately).

Weight: 3.5kg (typical model)

Capacity (total): 125ml

Consumption: 1ml of liquid gives approx., 200ml of vapour

Vapour flow pressure: At 5L/Min oxygen, 3cm w.g. "OFF" position.

Flow rate: Nominal 200ml - 15L/MIN

Compensation temperature range: 15° - 36°C

Connectors: Cagemount 23mm, BS22mm Uni-Fit or quick-release versions.

Filler: Keyed filler.

S/N: 049964

ADDITIONAL INFORMATION

Animal	N/A
SKU	600-001-cl
Fill Type	Key Fill
Fitting Type	Quick Release
anaesthesia_vap_anaesthetic	5% Isoflurane

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights