

MACH LED 130 WALL MOUNTED EXAMINATION WITH STERILIZABLE HANDLE - CLEARANCE

MACHLED130-W-SH-CL



PRODUCT DESCRIPTION

Mach LED 130 Wall Mounted Examination Light With Removable Sterilizable Handle

Producing a cool light the LED technology is much more effective than conventional light sources. The temperature increase in the surgeon's working area is almost non-existent. Features include a faceted multi-lens system, a colour rendering index of Ra=95 and superior handling design allowing the light to be positioned exactly to the wound field.

Specifications Common To All Models		
Technical Specifications	Mach LED 130 F	Mach LED 130
Light Intensity at 1 Meter Distance (lux)	60,000	50,000
Colour Rendering Index R(a) General	95	95
Light-Field Size (cm)	14 - 25 Focusable	Fixed Focus 17
Colour Temperature	4,500	4,500
Electronic Light Intensity Control On Lamp Head	50 - 100%	50 - 100%
Power Consumption	28 W	28 W
Number of LED's	19	19
Mains	230 V 50 Hz	230 V 50 Hz
Working Distance (cm)	70 - 140	70 - 140

Diameter of Light Head	33	33
Height Adjustment (cm)	123	123
Weight (Ceiling, Wall, Mobile)	15 Kg, 8 Kg, 20Kg	

S/N: 24-04613

ADDITIONAL INFORMATION

Animal	Mixed
SKU	MACHLED130-W-SH-cl
Autoclaveable Handle	Yes
Focusable	No
Head Type	Single
Light Source	LED
Lux @1 meter	60,000
Proceduure	Dental, Exam, Minor-Operation, Task

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights