

**WARMING TUBE FOR BURTONS EN-V5+ INFUSION PUMP**

EN-V5+WT

**PRODUCT DESCRIPTION****Add Infusion Warming To The EN-V5+ Pump With Ease**

Only available for the EN-V5+ Infusion Pump. The pump is compatible with the EN-V5WT warming lead and has a side holster built-in on the side of the unit to hold it. The lead operates up to 40°C to help keep your patient's temperature at healthy levels. The leads generous length allows for improved coverage and heat distribution across the IV giving set. The warming is only operational when running from mains power.

**ADDITIONAL INFORMATION**

<b>Animal</b>	Mixed
<b>SKU</b>	EN-V5+WT
<b>Catagory Sub Heading</b>	Only available for EN-V5+ Infusion Pump
<b>Flow Rate</b>	0.1 - 2000 ml/h

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights