

## SPITZ ALL IN ONE COMPLETE CAMERA SYSTEM

END-300-490-001



### PRODUCT DESCRIPTION

---

## A complete system for high-definition endoscopic, laparoscopic and arthroscopic surgery.

Universal and flexible "All in One Solution" for diagnostics and endoscopic therapy. The complete system can be easily transported, intuitively set up and used. The ultra-high resolution monitor and FullHD camera produce a high-quality image and the 60W equivalent light source provides good illumination.

#### System includes:

- Full HD Camera Head
- TV-Adapter, f1.4 to f3.5
- Light Guide cable with adapter
- Laptop with software
- LED Lightsource
- Suitcase

### ADDITIONAL INFORMATION

<b>Animal</b>	Companion animals
<b>SKU</b>	END-300-490-001
<b>Brand</b>	SPITZ

<b>Procedure Type</b>	Endoscopy, Laparoscopy, Arthroscopy, Cystoscopy, Tracheoscopy, Bronchoscopy, Cystoscopy, Otoscopy
<b>Endoscopy Device</b>	Camera System
<b>Animal Type</b>	Companion animals
<b>Image Resolution</b>	4K HD

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights