

SPITZ F1 CAMERA HEAD 1 HD

END-300-001-602

**PRODUCT DESCRIPTION****SPITZ F1 Camera Head 1 HD For use with the SPITZ F2 Pro Full HD Camera System**

Introducing the revolutionary SPITZ range, offering comprehensive solutions for veterinary diagnostics and procedures. The cutting-edge range encompasses diagnostic scopes, accessories, and everything required for rigid endoscopy procedures, including arthroscopy and laparoscopy.

ADDITIONAL INFORMATION

Animal	Companion animals, Equine and Large Animal
SKU	END-300-001-602
Brand	SPITZ
Procedure Type	Minimally invasive surgery, Endoscopy, Laparoscopy, Arthroscopy, Cystoscopy, Tracheoscopy, Bronchoscopy, Cystoscopy, Otoscopy
Endoscopy Device	Camera head
Animal Type	Companion animals, Equine and Large Animal
Image Resolution	Full HD

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights