

## SPITZ F1 CAMERA HEAD 1 HD

END-300-001-602



### PRODUCT DESCRIPTION

## SPITZ F1 Camera Head 1 HD For use with the SPITZ F2 Pro Full HD Camera System

Introducing the revolutionary SPITZ range, offering comprehensive solutions for veterinary diagnostics and procedures. The cutting-edge range encompasses diagnostic scopes, accessories, and everything required for rigid endoscopy procedures, including arthroscopy and laparoscopy.

### ADDITIONAL INFORMATION

<b>Animal</b>	Companion animals, Equine and Large Animal
<b>SKU</b>	END-300-001-602
<b>Brand</b>	SPITZ
<b>Procedure Type</b>	Minimally invasive surgery, Endoscopy, Laparoscopy, Arthroscopy, Cystoscopy, Tracheoscopy, Bronchoscopy, Cystoscopy, Otoscopy
<b>Endoscopy Device</b>	Camera head
<b>Animal Type</b>	Companion animals, Equine and Large Animal
<b>Image Resolution</b>	Full HD

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights