

SPITZ TV-ADAPTER F=25MM

END-300-001-814

**PRODUCT DESCRIPTION****SPiTZ Soakable C-mount video lense TV Coupler****Key features:**

- f=25mm manual lock soakable Full HD

ADDITIONAL INFORMATION

Animal	Companion animals, Equine and Large Animal
SKU	END-300-001-814
Brand	SPITZ
Procedure Type	Minimally invasive surgery, Endoscopy, Laparoscopy, Arthroscopy, Cystoscopy, Tracheoscopy, Bronchoscopy, Cystoscopy, Otoscopy
Endoscopy Device	C-Mount, video lens, TV Coupler
Animal Type	Companion animals, Equine and Large Animal
Image Resolution	Full HD

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights