

**SPITZ ZOOM-TV-ADAPTER F = 15.0MM - F = 36.0MM C-MOUNT**

END-250-015-036

**PRODUCT DESCRIPTION****SPiTZ Soakable C-Mount Zoom TV Coupler****Key Features:**

- f=15-36mm soakable Full HD

**ADDITIONAL INFORMATION**

<b>Animal</b>	Companion animals, Equine and Large Animal
<b>SKU</b>	END-250-015-036
<b>Brand</b>	SPiTZ
<b>Procedure Type</b>	Minimally invasive surgery, Endoscopy, Laparoscopy, Arthroscopy, Cystoscopy, Tracheoscopy, Bronchoscopy, Cystoscopy, Otoscopy
<b>Endoscopy Device</b>	C-Mount, video lens, TV Coupler
<b>Animal Type</b>	Companion animals, Equine and Large Animal
<b>Image Resolution</b>	Full HD

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights