

## SPITZ LUMEN LED 1 LIGHT SOURCE SYSTEM

END-300-900-180



### PRODUCT DESCRIPTION

---

## SPiTZ Lumen LED1 Light Source System

#### Key Features:

- Daylight intensity display
- 20.000h lamp lifetime
- Manual or automatic light intensity control
- Waterproof touch control panel
- Low heat conduction by improved LED technology

#### System includes:

- Main Unit
- BNC cable
- Mains cable
- Operating Manual

### ADDITIONAL INFORMATION

<b>Animal</b>	Companion animals, Equine and Large Animal
<b>SKU</b>	END-300-900-180
<b>Brand</b>	SPITZ

<b>Procedure Type</b>	Minimally invasive surgery, Endoscopy, Laparoscopy, Arthroscopy, Cystoscopy, Tracheoscopy, Bronchoscopy, Cystoscopy, Otoscopy
<b>Endoscopy Device</b>	Light source
<b>Animal Type</b>	Companion animals, Equine and Large Animal
<b>Image Resolution</b>	LED

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights