

## SPITZ ARTHRO P 1 IRRIGATION PUMP

END-225-906-000



### PRODUCT DESCRIPTION

## SPITZ ARTHRO P 1 Irrigation Pump For Arthroscopy

#### Key features:

- Hydrostatic pressure regulation Silicone Tubing Set included (END-225-910-001)
- Easy to operate, splash-proofed display
- Highly reliable tube guide and handling
- Skid-proof pedestals

#### Unit contains:

- Main Unit
- Tubing Set
- Mains Cable
- 2 Spare Fuses
- Operating Manual

### ADDITIONAL INFORMATION

|                         |  |
|-------------------------|--|
| <b>Animal</b>           | Companion animals, Equine and Large Animal |
| <b>SKU</b>              | END-225-906-000                            |
| <b>Brand</b>            | SPITZ                                      |
| <b>Procedure Type</b>   | Arthroscopy                                |
| <b>Endoscopy Device</b> | Arthroscopy device                         |
| <b>Animal Type</b>      | Companion animals, Equine and Large Animal |

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights