

**SPITZ AIRFLOW A 2 INSUFFLATOR 20L/MIN**

END-300-900-020

**PRODUCT DESCRIPTION**

**User friendly, highest safety standard for patients, pre-insfflation mode for use with Veress Cannulas, High gasflow rates**

**ADDITIONAL INFORMATION**

<b>Animal</b>	Companion animals, Equine and Large Animal
<b>SKU</b>	END-300-900-020
<b>Brand</b>	SPITZ
<b>Procedure Type</b>	Minimally invasive surgery, Laparoscopy
<b>Endoscopy Device</b>	Laparoscopy Device, Insufflator
<b>Animal Type</b>	Companion animals, Equine and Large Animal

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights