

## SPITZ LITHOTRIPTER L 100 UNIT

END-253-750-100



## PRODUCT DESCRIPTION

---

### SPiTZ Lithotripter L 100 Unit

#### Key Features:

- Energy transmission occur purely mechanically
- The smallest probe diameter is only 0.6mm, which ensures a larger lumen.
- Auxillary instruments can also be inserted.
- Very simple handling, all parts are easy to disassemble for cleaning purposes.
- Compatible to SPiTZ and other instruments for the therapy of ureteric, renal and bladder stones.
- Non-conductive hand piece to ensure safety for patients as well as users against current conduct.

#### System includes:

- Main Unit
- Mains cable
- Operating manual

## ADDITIONAL INFORMATION

<b>Animal</b>	Companion animals
<b>SKU</b>	END-253-750-100
<b>Brand</b>	SPiTZ

**Procedure Type**

Minimally invasive surgery

**Endoscopy Device**

Lithotriper

**Animal Type**

Companion animals

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights